

Anti-FANCA Antibody

Catalog Number: A03662-3

About FANCA

FANCA, also named as FAA, FACA and FANCH, is a DNA repair protein that may operate in a postreplication repair or a cell cycle checkpoint function. It's involved in interstrand DNA cross-link repair and in the maintenance of normal chromosome stability.

Overview

Product Name	Anti-FANCA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-FANCA Antibody catalog # A03662-3. Tested in WB, IHC, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O15360

Technical Details

Immunogen	E.coli-derived human FANCA recombinant protein (Position: D3-I294).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FANCA Antibody

For Research Use Only. Not for use in diagnostic procedures.