

## Anti-Cubilin CUBN Antibody

Catalog Number: A03670

### Overview


|                      |  |
|----------------------|--|
| Product Name         | Anti-Cubilin CUBN Antibody   |
| Reactive Species     | Human  |
| Description          | Boster Bio Anti-Cubilin CUBN Antibody catalog # A03670. Tested in WB,IHC applications. This antibody reacts with Human.                    |
| Application          | IHC, WB  |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | O60494   |

### Technical Details

|                     |                              |
|---------------------|------------------------------|
| Immunogen           | Recombinant protein          |
| Isotype             | IgG                          |
| Form                | Liquid                       |
| Concentration       | 1 mg/ml                      |
| Purification        | ProA affinity purified       |
| Suggested Dilutions | WB: 1:500<br>IHC: 1:50-1:200 |

## Anti-Cubilin CUBN Antibody (A03670) Images

---



Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-Cubilin antibody. Counter stained with hematoxylin.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cubilin CUBN Antibody

For Research Use Only. Not for use in diagnostic procedures.