

## Anti-ICAD (Q177) DFFA Antibody

Catalog Number: A03671-3

### Overview

Product Name	Anti-ICAD (Q177) DFFA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ICAD (Q177) DFFA Antibody catalog # A03671-3. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O00273

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 151-200 of Human ICAD.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

## Anti-ICAD (Q177) DFFA Antibody (A03671-3) Images




Figure 2. Immunohistochemistry validation of DFFA using Anti-ICAD (Q177) DFFA Antibody (A03671-3).

Immunohistochemistry (IHC) analyzes of ICAD (Q177) pAb in paraffin-embedded human breast carcinoma tissue at 1:100. For more protocol information of IHC




Figure 1. Western blotting validation for Anti-ICAD (Q177) DFFA Antibody A03671-3

Western blot (WB) analysis of ICAD (Q177) polyclonal antibody at 1:500 dilution

Lane1:C6 whole cell lysate(40ug)

Lane2:CT26 whole cell lysate(40ug)

Lane3:HCT116 whole cell lysate(40ug)

Lane4:COS-7 whole cell lysate(40ug)

Lane5:SGC7901 whole cell lysate(40ug)

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ICAD (Q177) DFFA Antibody