

## Anti-ORC5L/ORC5 Antibody Picoband® HRP Conjugated

Catalog Number: A03676-1-HRP

### About ORC5

Origin recognition complex subunit 5 is a protein that in humans is encoded by the ORC5 (ORC5L) gene. The origin recognition complex (ORC) is a highly conserved six subunit protein complex essential for the initiation of the DNA replication in eukaryotic cells. Studies in yeast demonstrated that ORC binds specifically to origins of replication and serves as a platform for the assembly of additional initiation factors such as Cdc6 and Mcm proteins. The protein encoded by this gene is a subunit of the ORC complex. Alternatively spliced transcript variants encoding distinct isoforms have been described.

### Overview

Product Name	Anti-ORC5L/ORC5 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Monkey, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O43913

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human ORC5, identical to the related mouse sequence.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ORC5L/ORC5 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.