

Anti-hnRNP U Antibody

Catalog Number: A03691-2


Overview

Product Name	Anti-hnRNP U Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-hnRNP U Antibody catalog # A03691-2. Tested in WB,ICC,IF,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q00839


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 ICC: 1:100-1:500 IHC: 1:50-1:200

Anti-hnRNP U Antibody (A03691-2) Images



Western blot analysis of HnRNP U on different cell lysates using anti-HnRNP U at 1/500 dilution. Positive control □ Lane 1: NIH-3T3 Lane 2: SH-SY-5Y Lane 3: K562



ICC staining HnRNP U in A431 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-hnRNP U Antibody

For Research Use Only. Not for use in diagnostic procedures.