

## Anti-Cytokeratin 9 KRT9 Antibody

Catalog Number: A03692-1

### Overview


Product Name	Anti-Cytokeratin 9 KRT9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cytokeratin 9 KRT9 Antibody catalog # A03692-1. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35527

### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200

## Anti-Cytokeratin 9 KRT9 Antibody (A03692-1) Images

---



Immunohistochemical analysis of paraffin-embedded mouse skin tissue using anti-Cytokeratin 9 antibody. Counter stained with hematoxylin.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytokeratin 9 KRT9 Antibody

For Research Use Only. Not for use in diagnostic procedures.