

Anti-Rab 18 Antibody

Catalog Number: A03719

About RAB18

CAND1 is also known as TIP120A, and TATA-binding protein-interacting protein 120A. The SCF complex consists of the invariable components Skp1, Cul1, and Rbx1 as well as a variable F-box protein, and functions as an E3 ubiquitin ligase. E3 ubiquitin ligases regulate various physiological processes. CAND1 binds to Cul1 and potentially regulates the SCF complex. CAND1 physically associates with Cul1 in the nucleus and this interaction is mediated by a central region of Cul1 distinct from its binding sites for Skp1 and Rbx1. CAND1 selectively binds to unneddylated CUL1 and is dissociated by CUL1 neddylation. CAND1 forms a ternary complex with CUL1 and ROC1.

Overview

Product Name	Anti-Rab 18 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Rab 18 Antibody catalog # A03719. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal 5A6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NP72

Technical Details

Immunogen	Synthesized peptide derived from the Internal region of human Rab 18. at AA range: 80-160
Predicted Reactive Species	Bovine, Canine, Equine, Guinea Pig, Pig, Rabbit
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this

kit.


If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.

Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:

WB, 1:500-1:2000

Anti-Rab 18 Antibody (A03719) Images



Western Blot (WB) analysis of MCF-7 cells using Rab 18 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rab 18 Antibody