

Anti-AChR Alpha 1 (E217) CHRNA1 Antibody

Catalog Number: A03723


Overview

Product Name	Anti-AChR Alpha 1 (E217) CHRNA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-AChR Alpha 1 (E217) CHRNA1 Antibody catalog # A03723. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P02708


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 181-230 of Human AChRalpha1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200


Anti-AChR Alpha 1 (E217) CHRNA1 Antibody (A03723) Images



Western blot (WB) analysis of AChRalpha1 (E217) pAb at 1:500 dilution Lane1:HepG2 whole cell lysate(40ug) Lane2:A549 whole cell lysate(40ug) Lane3:PC3 whole cell lysate(40ug) Lane4:CT26 whole cell lysate(40ug) Lane5:H9C2 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of AChRalpha1 (E217) pAb in paraffin-embedded human liver carcinoma tissue at 1:50, showing membrane and cytoplasmic staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.



Immunohistochemistry (IHC) analyzes of AChRalpha1 (E217) pAb in paraffin-embedded human kidney carcinoma tissue at 1:50, showing membrane and cytoplasmic staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AChR Alpha 1 (E217) CHRNA1 Antibody

For Research Use Only. Not for use in diagnostic procedures.