

Anti-GPR68 (R214) Antibody

Catalog Number: A03726-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-GPR68 (R214) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GPR68 (R214) Antibody catalog # A03726-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15743

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 183-232 of Human GPR68.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-GPR68 (R214) Antibody (A03726-1) Images



Western blot (WB) analysis of GPR68 (R214) pAb at 1:500 dilution Lane1:U-87MG whole cell lysate(40ug) Lane2:K562 whole cell lysate(40ug) Lane3:MCF-7 whole cell lysate(40ug) Lane4:The Lung tissue lysate of Mouse(40ug) Lane5:The Lung tissue lysate of Rat(40ug)

[Anti-GPR68 \(R214\) Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.