

Anti-ABCAC ABCA12 Antibody

Catalog Number: A03728

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-ABCAC ABCA12 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ABCAC ABCA12 Antibody catalog # A03728. Tested in IHC, IF applications. This antibody reacts with Human, Mouse.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86UK0

Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 2170-2250
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ABCAC antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:50-300 IF 1:50-200

Anti-ABCAC ABCA12 Antibody (A03728) Images



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).

[Anti-ABCAC ABCA12 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.