

Anti-AMPK Beta 1 PRKAB1 Antibody

Catalog Number: A03741-1


Overview

Product Name	Anti-AMPK Beta 1 PRKAB1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-AMPK Beta 1 PRKAB1 Antibody catalog # A03741-1. Tested in WB, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y478


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 FC: 1:50-1:100

Anti-AMPK Beta 1 PRKAB1 Antibody (A03741-1) Images



Western blot analysis of AMPK beta 1 on A431 cells lysates using anti-AMPK beta 1 antibody at 1/1,000 dilution.



Flow cytometric analysis of A431 cells with AMPK beta 1 antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AMPK Beta 1 PRKAB1 Antibody

For Research Use Only. Not for use in diagnostic procedures.