

Anti-BPTF/Falz Monoclonal Antibody

Catalog Number: A03754

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-BPTF/Falz Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-BPTF/Falz Monoclonal Antibody catalog # A03754. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 2F1
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q12830

Technical Details

Immunogen	Purified recombinant fragment of human BPTF expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000 This antibody reacts with human, rat, and mouse CAND1 tested by western blot and immunoprecipitation. The antibody

immunoprecipitates in vitro translated protein and protein from cell lysates (using HeLa and NIH-3T3, and others).
Coimmunoprecipitation of related proteins has not been tested. A 136.4 kDa band corresponding to human CAND1 is detected. Most cell lines expressing CAND1 can be used as a positive control. Researchers should determine optimal titers for other applications.

Anti-BPTF/Falz Monoclonal Antibody (A03754) Images



Western Blot analysis using BPTF Monoclonal Antibody against HEK293 (1) and BPTF (AA: 503-670)-hIgGfc transfected HEK293 (2) cell lysate.

Anti-BPTF/Falz Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.