

Anti-Internexin Alpha INA Antibody

Catalog Number: A03756


Overview

Product Name	Anti-Internexin Alpha INA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Internexin Alpha INA Antibody catalog # A03756. Tested in WB,ICC/IF,IHC,Flow Cytometry,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16352


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-Internexin Alpha INA Antibody (A03756) Images



Western blot analysis of alpha Internexin on different lysates using anti-alpha Internexin antibody at 1/1,000 dilution.
Positive control: Lane 1: Mouse brain tissue Lane 2: A549



ICC staining alpha Internexin in PC-12 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Internexin Alpha INA Antibody

For Research Use Only. Not for use in diagnostic procedures.