

## Anti-Integrin Beta 6/ITGB6 antibody

Catalog Number: A03800

### Overview

Product Name	Anti-Integrin Beta 6/ITGB6 antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Integrin Beta 6/ITGB6 antibody catalog # A03800. Tested in IHC, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P18564

### Technical Details

Immunogen	E.coli-derived human ITGB6 recombinant protein (Position: L141-S430)
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin Beta 6/ITGB6 antibody

For Research Use Only. Not for use in diagnostic procedures.