

Anti-Histamine N-methyltransferase HNMT Antibody

Catalog Number: A03807

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Histamine N-methyltransferase HNMT Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Histamine N-methyltransferase HNMT Antibody catalog # A03807. Tested in WB,IP applications. This antibody reacts with Human,Mouse.
Application	IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P50135


Technical Details

Immunogen	Recombinant fusion protein of human HNMT(NP_008826.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IP: 1:20-1:50

Anti-Histamine N-methyltransferase HNMT Antibody (A03807) Images



Western blot analysis of extracts of various cell lines, using HNMT antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunoprecipitation analysis of 200ug extracts of HT-29 cells using 1ug HNMT antibody. Western blot was performed from the immunoprecipitate using HNMT antibody at a dilution of 1:1000.

Anti-Histamine N-methyltransferase HNMT Antibody

For Research Use Only. Not for use in diagnostic procedures.