

## Anti-DDX6 Antibody

Catalog Number: A03826


### Overview

Product Name	Anti-DDX6 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DDX6 Antibody catalog # A03826. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P26196


### Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:100

## Anti-DDX6 Antibody (A03826) Images



Western blot analysis of DDX6 on different lysates using anti-DDX6 antibody at 1/500 dilution. Positive control: Lane 1: Human colon Lane 2: SH-SY-5Y Lane 3: SiHa Lane 4: K562



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-DDX6 antibody. Counter stained with hematoxylin.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-DDX6 Antibody

For Research Use Only. Not for use in diagnostic procedures.