

Anti-GABA B Receptor 2 Gabbr2 Antibody

Catalog Number: A03830-1


Overview

| | |
|----------------------|--|
| Product Name | Anti-GABA B Receptor 2 Gabbr2 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-GABA B Receptor 2 Gabbr2 Antibody catalog # A03830-1. Tested in WB,IHC,Flow Cytometry,ICC applications. This antibody reacts with Human,Mouse,Rat. |
| Application | Flow Cytometry, IHC, ICC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q80T41 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Recombinant protein |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | ProA affinity purified |
| Suggested Dilutions | WB: 1:500-1:1,000 IHC: 1:50-1:200 FC: 1:50-1:100 ICC: 1:50-1:100 |

Anti-GABA B Receptor 2 Gabbr2 Antibody (A03830-1) Images



Western blot analysis of GABA B Receptor 2 on mouse cerebellum tissue lysate using anti- GABA B Receptor 2 antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue using anti- GABA B Receptor 2 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GABA B Receptor 2 Gabbr2 Antibody

For Research Use Only. Not for use in diagnostic procedures.