

Anti-CD108 SEMA7A Antibody

Catalog Number: A03832


Overview

Product Name	Anti-CD108 SEMA7A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CD108 SEMA7A Antibody catalog # A03832. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75326

Technical Details

Immunogen	Synthesized peptide derived from CD108 . at AA range: 150-230
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:5000

Anti-CD108 SEMA7A Antibody (A03832) Images



Western blotting validation for Anti-CD108 SEMA7A Antibody A03832 Western Blot (WB) analysis of specific cells using CD108 polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD108 SEMA7A Antibody

For Research Use Only. Not for use in diagnostic procedures.