

Anti-USP22 Antibody

Catalog Number: A03887


Overview

| | |
|----------------------|--|
| Product Name | Anti-USP22 Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-USP22 Antibody catalog # A03887. Tested in WB,ICC/IF,Flow Cytometry,IP applications. This antibody reacts with Human. |
| Application | Flow Cytometry, IP, IF, ICC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q9UPT9 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Recombinant protein |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | ProA affinity purified |
| Suggested Dilutions | WB: 1:1,000-1:5,000 ICC: 1:50-1:200 FC: 1:50-1:100 IP: 1:10-1:50 |

Anti-USP22 Antibody (A03887) Images



Western blot analysis of USP22 on HeLa (1) and HepG2 (2) cell lysate using anti-USP22 antibody at 1/500 dilution.



ICC staining USP22 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-USP22 Antibody

For Research Use Only. Not for use in diagnostic procedures.