

## Anti-LAL LIPA Monoclonal Antibody

Catalog Number: A03895


### Overview

Product Name	Anti-LAL LIPA Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-LAL LIPA Monoclonal Antibody catalog # A03895. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 9G7F12; 7G6D7
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P38571

### Technical Details

Immunogen	Purified recombinant fragment of LAL expressed in E. Coli.
Isotype	IgG2a
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-LAL LIPA Monoclonal Antibody (A03895) Images



Western Blot analysis using LAL Monoclonal Antibody against LAL recombinant protein.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LAL LIPA Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.