

Anti-Histone acetyltransferase KAT7 KAT7 Antibody

Catalog Number: A03896-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Histone acetyltransferase KAT7 KAT7 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Histone acetyltransferase KAT7 KAT7 Antibody catalog # A03896-1. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95251


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 IHC: 1:50-1:100 RRM2B/p53R2 antibody is tested in Western Blot and ELISA and suitable for IP, IHC, and IF. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 41 kDa in size corresponding to RRM2B1 by western blotting in the appropriate cell lysate or extract.

Anti-Histone acetyltransferase KAT7 KAT7 Antibody (A03896-1) Images



Immunohistochemical analysis of paraffin-embedded mouse testes tissue using anti-KAT7 antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded rat testes tissue using anti-KAT7 antibody. Counter stained with hematoxylin.

Anti-Histone acetyltransferase KAT7 KAT7 Antibody

For Research Use Only. Not for use in diagnostic procedures.