

Anti-PKA R2 (Phospho-S99) PRKAR2A Antibody

Catalog Number: A03911S99


Overview

Product Name	Anti-PKA R2 (Phospho-S99) PRKAR2A Antibody
Reactive Species	Human, Mouse, Pig, Rat
Description	Boster Bio Anti-PKA R2 (Phospho-S99) PRKAR2A Antibody catalog # A03911S99. Tested in WB,ICC/IF,IHC,IP applications. This antibody reacts with Human,Mouse,Rat,Pig.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P13861


Technical Details

Immunogen	Synthetic phospho-peptide corresponding to residues surrounding Ser99 of human PKA R2
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 ICC: 1:50-1:200 IHC: 1:50-1:200

Anti-PKA R2 (Phospho-S99) PRKAR2A Antibody (A03911S99) Images



ICC staining phospho-PKA R2 (S99) in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



ICC staining phospho-PKA R2 (S99) in MCF-7 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PKA R2 (Phospho-S99) PRKAR2A Antibody

For Research Use Only. Not for use in diagnostic procedures.