

Anti-CD177 antigen CD177 Antibody

Catalog Number: A03958


Overview

Product Name	Anti-CD177 antigen CD177 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-CD177 antigen CD177 Antibody catalog # A03958. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8N6Q3


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human CD177. AA range:361-410
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-300 ELISA 1:20000 IF 1:50-200


Anti-CD177 antigen CD177 Antibody (A03958) Images



Western Blot analysis of K562 cells using CD177 Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Western blot analysis of lysate from K562 cells, using CD177 Antibody.



Immunohistochemical analysis of paraffin-embedded human colon, antibody was diluted at 1:100

1 Publications Citing This Product

1. PubMed ID: 10.1111/jpn.13664, Chitosan nanoparticles attenuate intestinal damage and inflammatory responses in LPS-challenged weaned piglets via prevention of IkappaB degradation

Visit bosterbio.com/anti-cd177-antibody-a03958-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD177 antigen CD177 Antibody

For Research Use Only. Not for use in diagnostic procedures.