

## Anti-p60 CAF1 CHAF1B Antibody

Catalog Number: A03959-1


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-p60 CAF1 CHAF1B Antibody   |
| Reactive Species     | Human   |
| Description          | Boster Bio Anti-p60 CAF1 CHAF1B Antibody catalog # A03959-1. Tested in WB,ICC,IF,IP,IHC,Flow Cytometry applications. This antibody reacts with Human. |
| Application          | Flow Cytometry, IP, IF, IHC, ICC, WB  |
| Clonality            | Polyclonal  |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2  |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.            |
| Host                 | Rabbit  |
| Uniprot ID           | Q13112  |


### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Recombinant protein corresponding to the C-terminus of Human p60 CAF1.                        |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | ProA affinity purified  |
| Suggested Dilutions | WB: 1:500-1:2,000<br>ICC/IF: 1:50-1:200<br>IHC: 1:50-1:200<br>IP: 1:10-1:50<br>FC: 1:50-1:100 |

## Anti-p60 CAF1 CHAF1B Antibody (A03959-1) Images



Western blot analysis of p60 CAF1 on different cell lysates using anti-p60 CAF1 antibody at 1/1,000 dilution. Positive control: Lane 1: SiHa Lane 2: K562 Lane 3: A431



ICC staining p60 CAF1 in 293T cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p60 CAF1 CHAF1B Antibody

For Research Use Only. Not for use in diagnostic procedures.