

Anti-Protein BTG1 BTG1 Antibody

Catalog Number: A03967


Overview

Product Name	Anti-Protein BTG1 BTG1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Protein BTG1 BTG1 Antibody catalog # A03967. Tested in WB,ICC/IF applications. This antibody reacts with Human.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P62324

Technical Details

Immunogen	Recombinant fusion protein of human BTG1(NP_001722.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:100

Anti-Protein BTG1 BTG1 Antibody (A03967) Images



Immunofluorescence analysis of U2OS cells using BTG1 antibody .

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein BTG1 BTG1 Antibody

For Research Use Only. Not for use in diagnostic procedures.