

Anti-Gasdermin-B GSDMB Antibody

Catalog Number: A03968


Overview

Product Name	Anti-Gasdermin-B GSDMB Antibody
Reactive Species	Human
Description	Boster Bio Anti-Gasdermin-B GSDMB Antibody catalog # A03968. Tested in WB,ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8TAX9

Technical Details

Immunogen	Recombinant full length protein of human GSDML.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ELISA: 1:5000-1:10000

Anti-Gasdermin-B GSDMB Antibody (A03968) Images



Western blot analysis of GSDML expression in HeLa (A), HT29 (B), HL-60 (C), MDA-MB435 (D) whole cell lysates.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Gasdermin-B GSDMB Antibody

For Research Use Only. Not for use in diagnostic procedures.