

Anti-UBE2D2/Ubch5B Antibody

Catalog Number: A03971-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-UBE2D2/Ubch5B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-UBE2D2/Ubch5B Antibody catalog # A03971-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P62837

Technical Details


Immunogen	The antiserum was produced against synthesized peptide derived from human UBE2D2. AA range:98-147
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000

	ELISA 1:10000
--	---------------

Anti-UBE2D2/Ubch5B Antibody (A03971-1) Images



Western Blot analysis of various cells using UBE2D2 Polyclonal Antibody diluted at 1:2000. Secondary antibody was diluted at 1:20000



Western blot analysis of the lysates from K562 cells using UBE2D2 antibody.

Anti-UBE2D2/Ubch5B Antibody

For Research Use Only. Not for use in diagnostic procedures.