

Anti-PPP1CB Antibody

Catalog Number: A04022-2


Overview

Product Name	Anti-PPP1CB Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PPP1CB Antibody catalog # A04022-2. Tested in WB,IP,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P62140


Technical Details

Immunogen	Recombinant protein corresponding to human PPP1CB aa 100-327.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified.
Suggested Dilutions	WB: 1:1,000-1:10,000 IP: 1:50 IHC: 1:50-1:200

Anti-PPP1CB Antibody (A04022-2) Images



Western blot analysis of PPP1CB on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody was used at a 1:500 dilution in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature. Positive control: Lane 1: Mouse skeletal muscle tissue lysate Lane 2: Rat brain tissue lysate Lane 3: A431 cell lysate Lane 4: SH-SY-5Y cell lysate Lane 5: SiHa cell lysate Lane 6: 293 cell lysate Lane 7: Hela cell lysate Lane 8: PC-3M cell lysate Lane 9: HepG2 cell lysate Lane 10: A549 cell lysate



Immunohistochemical analysis of paraffin-embedded rat brain tissue using anti-PPP1CB antibody. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0) for 20 mins. The tissues were blocked in 5% BSA for 30 minutes at room temperature, washed with ddH₂O and PBS, and then probed with ET7109-40 at 1/200 dilution, for 30 minutes at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. Counter stained with hematoxylin and mounted with DPX.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PPP1CB Antibody

For Research Use Only. Not for use in diagnostic procedures.