

Anti-BOP1 Antibody Picoband® HRP Conjugated

Catalog Number: A04026-2-HRP

About BOP1

Enables RNA binding activity. Involved in regulation of cell cycle; regulation of signal transduction by p53 class mediator; and ribosomal large subunit biogenesis. Located in chromosome; nucleolus; and nucleoplasm. Part of PeBoW complex.

Overview

Product Name	Anti-BOP1 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q14137

Technical Details

Immunogen	E.coli-derived human BOP1 recombinant protein (Position: R157-N701). Human BOP1 shares 89.9% and 89.7% amino acid (aa) sequence identity with mouse and rat BOP1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BOP1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.