

Anti-KIR2DS4 Antibody

Catalog Number: A04052


Overview

Product Name	Anti-KIR2DS4 Antibody
Reactive Species	Human
Description	Boster Bio Anti-KIR2DS4 Antibody catalog # A04052. Tested in WB,IF applications. This antibody reacts with Human.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P43632

Technical Details

Immunogen	Recombinant full length Human KIR2DS4.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IF: 1:50-1:200

Anti-KIR2DS4 Antibody (A04052) Images



Western blot analysis of extract of various cells, using KIR2DS4 antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KIR2DS4 Antibody

For Research Use Only. Not for use in diagnostic procedures.