

## Anti-ARH/LDLRAP1 Antibody

Catalog Number: A04054-2

### About LDLRAP1

The protein encoded by this gene is a cytosolic protein which contains a phosphotyrosine binding (PTD) domain. The PTD domain has been found to interact with the cytoplasmic tail of the LDL receptor. Mutations in this gene lead to LDL receptor malfunction and cause the disorder autosomal recessive hypercholesterolaemia.

### Overview

Product Name	Anti-ARH/LDLRAP1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ARH/LDLRAP1 Antibody catalog # A04054-2. Tested in WB, IHC, IP, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5SW96

### Technical Details

Immunogen	E.coli-derived human ARH/LDLRAP1 recombinant protein (Position: A3-F306).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ARH/LDLRAP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.