

Anti-Secretogranin-1 CHGB Antibody

Catalog Number: A04073


Overview

Product Name	Anti-Secretogranin-1 CHGB Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Secretogranin-1 CHGB Antibody catalog # A04073. Tested in WB applications. This antibody reacts with Human,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P05060

Technical Details

Immunogen	Recombinant fusion protein of human CHGB(NP_001810.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:100-1:500

Anti-Secretogranin-1 CHGB Antibody (A04073) Images



Western blot analysis of extracts of various cell lines, using CHGB antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Secretogranin-1 CHGB Antibody

For Research Use Only. Not for use in diagnostic procedures.