

## Anti-FRK Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A04075-1-Fluoro488

### About FRK

Fyn-related kinase (Frk) is a tyrosine kinase of the Src family. Also known as Rak, Gtk, and Iyk, Frk is predominantly expressed in epithelial tissues. Frk protein is 49% and 47% identical to Fyn and Csk, respectively. Src family kinases regulate an array of cellular processes, including growth factor signaling, cytoskeleton dynamics, and cell proliferation.

### Overview

Product Name	Anti-FRK Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P42685

### Technical Details

Immunogen	E.coli-derived human FRK recombinant protein (Position: Q39-A339).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FRK Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.