

## Anti-ARL3 Antibody Picoband® HRP Conjugated

Catalog Number: A04080-1-HRP

### About ARL3

ADP-ribosylation factor-like protein 3 is a protein that in humans is encoded by the ARL3 gene. Enables GDP binding activity; GTP binding activity; and microtubule binding activity. Involved in several processes, including cilium assembly; protein localization to cilium; and small GTPase mediated signal transduction. Acts upstream of or within post-Golgi vesicle-mediated transport. Located in several cellular components, including microtubule cytoskeleton; midbody; and photoreceptor connecting cilium. Implicated in Joubert syndrome and retinitis pigmentosa 83.

### Overview

Product Name	Anti-ARL3 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P36405

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human ARL3, identical to the related mouse and rat sequences.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ARL3 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.