

## Anti-CD319/SLAMF7 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A04083-2-Fluoro647

### About SLAMF7

Enables identical protein binding activity. Predicted to be involved in T cell activation and immune response. Predicted to act upstream of or within regulation of natural killer cell activation. Located in endoplasmic reticulum.

### Overview

Product Name	Anti-CD319/SLAMF7 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NQ25

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human CD319/SLAMF7.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-CD319/SLAMF7 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.