

## Anti-ACF1/BAZ1A Antibody Picoband® Biotin Conjugated

Catalog Number: A04090-1-Biotin

### About BAZ1A

Bromodomain adjacent to zinc finger domain protein 1A is a protein that in humans is encoded by the BAZ1A gene. The BAZ1A gene encodes the accessory subunit of the ATP-dependent chromatin assembly factor (ACF), a member of the ISWI ('imitation switch') family of chromatin remodeling complexes.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-ACF1/BAZ1A Antibody Picoband® Biotin Conjugated   |
| Reactive Species     | Human, Monkey, Mouse, Rat  |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.                             |
| Host                 | Rabbit   |
| Uniprot ID           | Q9NRL2   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human ACF1/BAZ1A recombinant protein (Position: E1033-H1536).                       |
| Cross Reactivity    | No cross-reactivity with other proteins.   |
| Isotype             | Rabbit IgG   |
| Form                | Liquid   |
| Concentration       | 0.5 mg/mL  |
| Purification        | Immunogen affinity purified.   |
| Conjugate           | Biotin   |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ACF1/BAZ1A Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.