

## Anti-IGFBP6 Antibody

Catalog Number: A04101

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



## Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                     |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-IGFBP6 Antibody                                                                                                                                                                                                                                                                                                                                                                                                |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                   |
| Description          | Boster Bio Anti-IGFBP6 Antibody catalog # A04101. Tested in ELISA, WB applications. This antibody reacts with Human.                                                                                                                                                                                                                                                                                                |
| Application          | ELISA, WB                                                                                                                                                                                                                                                                                                                                                                                                           |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                          |
| Formulation          | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                          |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                              |
| Uniprot ID           | P24592                                                                                                                                                                                                                                                                                                                                                                                                              |

## Technical Details

|               |                                                                                                                           |
|---------------|---------------------------------------------------------------------------------------------------------------------------|
| Immunogen     | The antiserum was produced against synthesized peptide derived from the Internal region of human IGFBP6. AA range:101-150 |
| Isotype       | IgG                                                                                                                       |
| Form          | Liquid                                                                                                                    |
| Concentration | 1 mg/ml                                                                                                                   |
| Purification  | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.     |

Suggested Dilutions

WB 1:500-1:2000  
ELISA 1:20000

Affinity purified Bach-2 antibody has been tested for use in ELISA, Immunoprecipitation, and western blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 91-140 kDa in size corresponding to Bach2 protein by western blotting in the appropriate cell lysate or extract.

## Anti-IGFBP6 Antibody (A04101) Images

---



Western blotting validation for Anti-IGFBP6 Antibody A04101  
Western Blot (WB) analysis of K562 cells using IGFBP6 polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

### Anti-IGFBP6 Antibody

For Research Use Only. Not for use in diagnostic procedures.