

## Anti-SFRP5 Antibody

Catalog Number: A04109

### About SFRP5

Secreted frizzled-related protein 5 (SFRP5) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. SFRP5 and SFRP1 may be involved in determining the polarity of photoreceptor cells in the retina. SFRP5 is highly expressed in the retinal pigment epithelium, and moderately expressed in the pancreas.

### Overview

Product Name	Anti-SFRP5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SFRP5 Antibody catalog # A04109. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt $-20^{\circ}\text{C}$ as supplied. 6 months $2$ to $8^{\circ}\text{C}$ after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5T4F7

### Technical Details

Immunogen	E.coli-derived human SFRP5 recombinant protein (Position: A3-A314).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SFRP5 Antibody

For Research Use Only. Not for use in diagnostic procedures.