

Anti-MST3/STK24 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A04112-3-Fluoro550

About STK24

Serine/threonine-protein kinase 24 is an enzyme that in humans is encoded by the STK24 gene located in the chromosome 13, band q32.2.

Overview

Product Name	Anti-MST3/STK24 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y6E0

Technical Details

Immunogen	E.coli-derived human MST3/STK24 recombinant protein (Position: M1-Y433).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MST3/STK24 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.