

Anti-LGP2/DHX58 Antibody Picoband® Biotin Conjugated

Catalog Number: A04121-2-Biotin

About DHX58

Enables several functions, including ATP hydrolysis activity; RNA binding activity; and zinc ion binding activity. Involved in negative regulation of type I interferon production and regulation of innate immune response. Predicted to be active in cytoplasm.

Overview

| | |
|----------------------|--|
| Product Name | Anti-LGP2/DHX58 Antibody Picoband® Biotin Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | Q96C10 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human LGP2/DHX58 recombinant protein (Position: A38-S658). Human LGP2/DHX58 shares 78.4% amino acid (aa) sequence identity with mouse LGP2/DHX58. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Biotin |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LGP2/DHX58 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.