

Anti-LGP2/DHX58 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A04121-2-Fluoro488

About DHX58

Enables several functions, including ATP hydrolysis activity; RNA binding activity; and zinc ion binding activity. Involved in negative regulation of type I interferon production and regulation of innate immune response. Predicted to be active in cytoplasm.

Overview

Product Name	Anti-LGP2/DHX58 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q96C10

Technical Details

Immunogen	E.coli-derived human LGP2/DHX58 recombinant protein (Position: A38-S658). Human LGP2/DHX58 shares 78.4% amino acid (aa) sequence identity with mouse LGP2/DHX58.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-LGP2/DHX58 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.