

Anti-LGP2/DHX58 Antibody Picoband® HRP Conjugated

Catalog Number: A04121-2-HRP

About DHX58

Enables several functions, including ATP hydrolysis activity; RNA binding activity; and zinc ion binding activity. Involved in negative regulation of type I interferon production and regulation of innate immune response. Predicted to be active in cytoplasm.

Overview

Product Name	Anti-LGP2/DHX58 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q96C10

Technical Details

Immunogen	E.coli-derived human LGP2/DHX58 recombinant protein (Position: A38-S658). Human LGP2/DHX58 shares 78.4% amino acid (aa) sequence identity with mouse LGP2/DHX58.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LGP2/DHX58 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.