

## Anti-HSD17B3 Antibody Picoband® Cy3 Conjugated

Catalog Number: A04128-3-Cy3

### About HSD17B3

17beta-Hydroxysteroid dehydrogenase 3 (17beta-HSD3) is an enzyme that in humans is encoded by the HSD17B3 gene and is involved in androgen steroidogenesis. This isoform of 17 beta-hydroxysteroid dehydrogenase is expressed predominantly in the testis and catalyzes the conversion of androstenedione to testosterone. It preferentially uses NADP as cofactor. Deficiency can result in male pseudohermaphroditism with gynecomastia.

### Overview

Product Name	Anti-HSD17B3 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P37058

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human HSD17B3, which shares 68.4% amino acid (aa) sequence identity with rat HSD17B3.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSD17B3 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.