

Anti-p23 PTGES3 Antibody

Catalog Number: A04136-3


Overview

Product Name	Anti-p23 PTGES3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-p23 PTGES3 Antibody catalog # A04136-3. Tested in WB,ICC,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q15185


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC: 1:50-1:200 IHC: 1:100-1:500 IP: 1:10-1:50

Anti-p23 PTGES3 Antibody (A04136-3) Images



Western blot analysis of p23 on different lysates using anti-p23 antibody at 1/500 dilution. Positive control: Lane 1: Mouse brain tissue Lane 2: SK-Br-3 Lane 3: Rat lung tissue



ICC staining p23 in A431 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p23 PTGES3 Antibody

For Research Use Only. Not for use in diagnostic procedures.