

Anti-Nuclear factor 1 X-type NFIX Antibody

Catalog Number: A04138


Overview

Product Name	Anti-Nuclear factor 1 X-type NFIX Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Nuclear factor 1 X-type NFIX Antibody catalog # A04138. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14938

Technical Details

Immunogen	Recombinant fusion protein of human NFIX(NP_002492.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Nuclear factor 1 X-type NFIX Antibody (A04138) Images



Western blot analysis of extracts of various cell lines, using NFIX antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nuclear factor 1 X-type NFIX Antibody

For Research Use Only. Not for use in diagnostic procedures.