

## Anti-ARL2 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A04155-1-Fluoro594

### About ARL2

ADP-ribosylation factor-like protein 2 is a protein that in humans is encoded by the ARL2 gene. This gene encodes a small GTP-binding protein of the RAS superfamily which functions as an ADP-ribosylation factor (ARF). The encoded protein is one of a functionally distinct group of ARF-like genes.

### Overview

Product Name	Anti-ARL2 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P36404

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human ARL2, which shares 90.9% amino acid (aa) sequence identity with mouse and rat ARL2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-ARL2 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.