

Anti-FKBP51 Antibody

Catalog Number: A04182-1


Overview

Product Name	Anti-FKBP51 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-FKBP51 Antibody catalog # A04182-1. Tested in WB,IHC,IP,Flow Cytometry applications. This antibody reacts with Human,Rat.
Application	Flow Cytometry, IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13451


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-FKBP51 Antibody (A04182-1) Images



Western blot analysis of FKBP51 on rat kidney tissue and Daudi cell lysate using anti-FKBP51 antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-FKBP51 antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FKBP51 Antibody

For Research Use Only. Not for use in diagnostic procedures.