

Anti-Myb-binding protein 1A MYBBP1A Antibody

Catalog Number: A04187


Overview

Product Name	Anti-Myb-binding protein 1A MYBBP1A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Myb-binding protein 1A MYBBP1A Antibody catalog # A04187. Tested in WB,IHC,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BQG0


Technical Details

Immunogen	Recombinant fusion protein of human MYBBP1A(NP_055335.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 ICC/IF: 1:50-1:200


Anti-Myb-binding protein 1A MYBBP1A Antibody (A04187) Images




Western blot analysis of extracts of various cell lines, using MYBBP1A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Confocal immunofluorescence analysis of U2OS cells using MYBBP1A Polyclonal Antibody at dilution of 1:200. Blue: DAPI for nuclear staining.



Immunohistochemistry of paraffin-embedded rat lung using MYBBP1A Rabbit pAb at dilution of 1:100. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded rat spleen using MYBBP1A Rabbit pAb at dilution of 1:100. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Myb-binding protein 1A MYBBP1A Antibody

For Research Use Only. Not for use in diagnostic procedures.