

Anti-Myb-binding protein 1A MYBBP1A Antibody

Catalog Number: A04187-1


Overview

Product Name	Anti-Myb-binding protein 1A MYBBP1A Antibody
Reactive Species	Human
Description	Boster Bio Anti-Myb-binding protein 1A MYBBP1A Antibody catalog # A04187-1. Tested in WB,Flow Cytometry,IHC,ICC,IF applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BQG0


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 FC: 1:50-1:10 IHC: 1:50 ICC: 1:10-1:50

Anti-Myb-binding protein 1A MYBBP1A Antibody (A04187-1) Images



Western blot analysis of MYBBP1A on A549 cell using anti-MYBBP1A antibody at 1/500 dilution.



Flow cytometric analysis of 293T cells with MYBBP1A antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti-rabbit IgG was used as the secondary antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Myb-binding protein 1A MYBBP1A Antibody

For Research Use Only. Not for use in diagnostic procedures.