

Anti-Laminnin Alpha 3 LAMA3 Antibody

Catalog Number: A04194

Overview

Product Name	Anti-Laminnin Alpha 3 LAMA3 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Laminnin Alpha 3 LAMA3 Antibody catalog # A04194. Tested in WB,ICC/IF applications. This antibody reacts with Human.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16787


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:100-1:500

Anti-Laminnin Alpha 3 LAMA3 Antibody (A04194) Images



Western blot analysis of Laminnin 5 alpha 3 on A431 cells lysates using anti-Laminnin 5 alpha 3 antibody at 1/1,000 dilution.



ICC staining Laminnin 5 alpha 3 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Laminnin Alpha 3 LAMA3 Antibody

For Research Use Only. Not for use in diagnostic procedures.